



LG4

LG4 – a versatile
radiofrequency
generator for pain
treatment




C3

C3 – a dedicated
cryoneurolysis system

Fire & Ice

inomed is the only company to offer both
radiofrequency & cryoanalgesia systems

inomed 



Humanity has been using cold to relieve pain for millennia



1977


Lloyd coined the term cryoanalgesia in a paper published in "The Lancet" claiming this technique was superior because it is not followed by neuritis or neuroma



2013

inomed introduces the first C3 CryoSystem

What is cryoanalgesia and how does it work?




4

Cryoanalgesia is a specialised technique for providing long-term pain relief where peripheral nerves are partially destroyed by extreme cold in order to create a conduction block and stop the pain signal.

Although cryoanalgesia is an old-established technique, it is gaining a resurgence as an interventional pain technique as the lesion size is comparatively much larger than a radiofrequency heat lesion. It is also reversible and hence a repeatable technique.


The ice crystals created during the procedure disrupt the axon causing Wallerian degeneration distal and slightly proximal to the lesion. The ice crystals only affect the axon and myelin sheath.



Cryoanalgesia is ideally suited to treat large myelinated sensory or mixed nerves that are responsible for peripheral pain

The cells regenerate at a rate of 1 mm per day and the analgesic effect is temporary and dependent on the quality of the freezing. Although the analgesic effect continues after the axon has regenerated, the duration of pain relief felt by the patient is comparable to that of a radiofrequency procedure.


The long term analgesic affect is due to the fact that ice crystals, which destroy the axon, also create vascular damage to the vasa nervorum, which produces a severe endoneural oedema.




1. Ice ball formation

An ice ball is formed around the targeted nerve, freezing induces either a neuropraxic or axonotmesis injury depending on the temperature attained.


Nerve vessels






2. Degeneration

When temperatures of more than -20°C are reached, Wallerian degeneration of the axon ensues.



3. Regeneration


The axon regenerates at a rate of 1 mm per day.




4. Reinnervation

The axon is fully regenerated. The analgesic effect is maintained even when the nerves have fully regenerated.

Structure of a nerve





James Prescott Joule



William Thomson


Two British physicists who developed the thermodynamic principles that the cryoprobe utilises in order to create an ice ball

Physics behind a cryoanalgesia cryoprobe

A cryoanalgesia cryoprobe works by using the Joule-Thomson effect: a compressed gas, either carbon dioxide (CO_2) or nitrous oxide (N_2O), is passed through a narrow tube and then suddenly allowed to expand at the tip of the outer sealed probe.

The pressure drop from 650–800 psi to 80–100 psi results in a rapid decrease in temperature and a cooling of the probe tip. Absorption of heat from surrounding tissues accompanies expansion of any gas, according to the principles of general gas law; this is the adiabatic principle of gas cooling extraction or the Joule-Thomson effect.

The tissue around the probe rapidly cools and an ice ball is formed. The ice ball varies in size depending on probe size, freeze time, tissue permeability to water and the presence of any vascular structures, which may cause a heat sink. Precise levels of gas flow through the cryoprobe are essential for maximum efficiency. Inadequate gas flow does not freeze tissue, excessive gas flow results in freezing down the stem of the cryoprobe.



Is it possible to cause permanent damage to the peripheral nerve?

The Sunderland classification of peripheral nerve injury, which was developed from the earlier Seddon classification, outlines the grading of nerve injuries. The application of extreme cold to a nerve root creates a 2nd degree injury or axonotmesis.

The C3 CryoSystem utilises carbon dioxide (CO₂), which has a cryo boiling point of -79 °C or nitrous oxide (N₂O), which has a cryo boiling point of -88 °C. This means that the C3 CryoProbe can never create a non-reversible cryolesion (Sunderland grade 3 and above), as the temperature at the probe tip cannot be lower than the cryogenic boiling point of the gas used to create the ice ball.

Reversible

1st degree

Neuropraxia – interruption of conduction, short recovery time

+10 °C to -20 °C

2nd degree

Axonotmesis – loss of continuity of the axon; Wallerian degeneration; preservation of endo-, peri- and epineurium

-20 °C to -100 °C
C3 CryoSystem

Non-reversible

3rd / 4th
degree

Neurotmesis – loss of continuity, and colder some loss of continuity of epineurium and perineurium


-140 °C

5th degree



Transection (severe neurotmesis) – gross loss of continuity



CRF isotherm showing temperature dissipation



Comparison of CRF heat lesion and cryo ice ball lesion

	Thermoprobe		Cryoprobe	
Probe diameter	0.7 mm Ø	0.9 mm Ø	1.3 mm Ø	2.1 mm Ø
				
Lesion diameter	2.4 mm Ø	3.1 mm Ø	3.5 mm Ø	5.5 mm Ø

What is the difference between RF lesion and cryolesion?


Continuous radiofrequency (CRF) lesion is achieved by prolonged application of continuous electrical current which heats the tissues to neurodestructive temperatures of 60–80 °C. This causes coagulative necrosis of both cellular and acellular structures. Histologically, both axonal degeneration and collagen fibre destruction of the endoneurium, perineurium and epineurium structures occurs. CRF produces a 3rd or 4th degree nerve injury with Wallerian degeneration and associated potential for neuroma and neuritis. CRF is a highly versatile technique, that creates relatively small lesions, which require a high degree of placement accuracy of the cannula and subsequent electrode placement.

Pulsed radiofrequency (PRF) induces a temporary electromagnetic field, which creates cellular changes that favourably alter the transmission of pain signals, although the true mechanism of how this is achieved remains unclear. The popular view is that PRF alters gene expression, neuronal membrane function and cytokine regulation. PRF does not fit into the Sunderland classification as it is not known how the nerve is altered to suppress the pain signal. PRF creates no tissue destruction and as such it does not have any risk of neuroma formation. The pain relief from PRF is of a shorter duration when compared to CRF or cryoanalgesia. In addition, the way in which PRF works is not sufficiently understood, which has led to significant differences in its application, calling into question the effectiveness of the procedure.

Cryoanalgesia utilises extreme cold to create a reversible pain conduction block. This is achieved by the formation of an ice ball around the target nerve. The ice crystals create an axonotmesis injury, which creates a reversible lesion that can be easily seen via ultrasound. Cryo ice ball lesions are comparatively much larger than CRF lesions and, as in the case of PRF, can also be used on mixed sensorimotor peripheral nerves, since cryoneurolysis does not cause any long-term motor deficits. This has led to cryoanalgesia being utilised as a treatment for spasticity as clinicians are able to treat pure motor nerves, pure sensory nerves and mixed sensorimotor nerves.

inomed is the only company that can offer radiofrequency and cryoanalgesia interventional pain treatment solutions





Unique benefits from cryoanalgesia treatment

- Direct visualisation of lesion under ultrasound
→ Can be performed in clinical setting
- In contrast to surgical, chemical or heat ablation, cryoneurolysis does not disrupt the epineurium, endoneurium or perineurium
→ No risk of procedural neuritis
- Reversible → Procedure can be safely repeated
- Low risk treatment option
- Wide range of clinical applications
- Large elliptical lesion compared to radiofrequency
- Immediate pain relief
- Pain relief from 6 months to 2 years
- Effective treatment option for spasticity
- Can be used preoperatively to reduce postoperative pain and opioids in arthroplasty

Clinical applications for cryoanalgesia

Craniofacial pain

- Supraorbital neuralgia
- Infraorbital neuralgia
- Mandibular neuralgia
- Mental neuralgia
- Auriculotemporal neuralgia
- Trigeminal neuralgia
- Posterior auricular neuralgia
- Glossopharyngeal neuralgia

Abdominal and pelvic

- Iliohypogastric neuralgia
- Ilioinguinal neuralgia
- Subgastric neuralgia
- Genitofemoral
- Sacral neuralgia
- Pudendal neuralgia
- Post-thoracotomy neuromas
- Persistent pain post rib fracture
- Thoracic post-herpetic neuralgia


Low back and lower extremity

- Facet joint pain
- Interspinous ligament
- Pseudosciatica
- Superior gluteal nerve
- Sacroiliac joint pain
- Cluneal neuralgia
- Lower extremity pain
- Obturator nerve
- Deep peroneal nerve
- Medial & lateral calcaneal nerves
- Peripheral neuropathy
- Phantom limb pain

Upper extremity pain

- Suprascapular nerve
- Superficial radial nerve

Indications Reference – Cryoanalgesia in interventional pain management. Pain Physician 2003 Jul;6(3):345-60
- Dr A. Trescot MD, DABBIPP, FIPP, CIPS – Past President of American Society of Interventional Pain Physicians.



Features of the C3 CryoSystem

The C3 CryoSystem generator is a compact single-channel device that features an easy-to-follow, intuitive graphic display. The C3 has been designed to be used with carbon dioxide (CO₂) and nitrous oxide (N₂O) from a non-siphoned gas cylinder.

The C3 has a list of key features that enhance the end users experience when delivering a cryoneurolysis treatment.

1

Automatic control and regulation of gas flow and pressure

2

Sensory and motor stimulation

3

Six soft key function buttons

4

Central rotary control knob

5

Neutral electrode connection

6

CryoProbe connection port

Additional features include


- Probe recognition
- USB data flash drive
- Barcode reader for patient identification
- Creation of treatment reports
- Acoustic signal to indicate freezing and defrosting cycle

Features of the C3 CryoSystem: C3 CryoProbes

The C3 CryoProbes have been specifically designed to deliver a targeted cryogenic ice ball to peripheral nerves for the purpose of providing prolonged pain relief.

The probes feature a coupling plug, which incorporates various types of connections to the C3 CryoSystem.



The C3 CryoProbes also have an integrated memory that contains the information about the probe type, the required probe pressure and usage information.



Features of the C3 CryoSystem: Foot switch

The C3 CryoSystem has an optional foot switch that enables the CryoSystem to be operated independently by the clinician. The foot switch is operated by short and long presses on the three switches.

The foot switch functions are displayed at the base of the C3 CryoSystem display screen.




In order to generate an ice ball, the C3 CryoProbes use the Joule-Thomson effect of gas expansion. Cryogenic gas is forced at high pressure to the tip of the probe, where it passes through a Joule-Thomson nozzle, allowing the gas to expand. The resultant pressure drop causes the cryogenic gas to cool rapidly, which freezes the water content in the adjacent tissue.

Rapid defrost is achieved by equalising the pressure on either side of the Joule-Thomson nozzle. Once the defrost is completed, the cryoprobe can be safely removed without any fear of tissue adhesion to the cryoprobe.

The C3 CryoProbes have a built-in nerve stimulator for precise targeting of the nerve. They are monopolar and require a neutral surface electrode positioned on the patient's body to complete the electrical circuit.

The C3 CryoProbes are reusable and available in two sizes: 1.3 mm blunt (15 gauge) and 2.1 mm sharp (12 gauge).





Cannulas

Art. No.	Name	Description
217604	Disposable cannula 15Gx45 with mandrin	Suitable for C3 CryoProbe 1.3 mm (217100) With depth marking for orientation Overall length 74.6 mm, working length 45 mm Outside diameter 1.8 mm (15G) Single use product, EO sterilised
217606	Disposable cannula 15Gx65 with mandrin	Suitable for C3 CryoProbe 1.3 mm (217100) With depth marking for orientation Overall length 94.6 mm, working length 65 mm Outside diameter 1.8 mm (15G) Single use product, EO sterilised
217609	Disposable cannula 15Gx95 with mandrin	Suitable for C3 CryoProbe 1.3 mm (217100) With depth marking for orientation Overall length 124.6 mm, working length 95 mm Outside diameter 1.8 mm (15G) Single use product, EO sterilised
217660	Disposable cannula 12Gx62 with mandrin	Suitable for C3 CryoProbe 2.1 mm (217200) With depth marking for orientation Overall length 96.3 mm, working length 62.5 mm Outside diameter 2.6 mm (12G) Single use product, EO sterilised
217680	Disposable cannula 12Gx92 with mandrin	Suitable for C3 CryoProbe 2.1 mm (217200) With depth marking for orientation Overall length 126.3 mm, working length 92.5 mm Outside diameter 2.6 mm (12G) Single use product, EO sterilised



CryoProbes



Art. No.	Name	Description	
217100	C3 CryoProbe 1.3 mm blunt for C3 CryoSystem	Diameter 1.3 mm, blunt tip Stem length 133 mm Handle length 140 mm Hose length 2 m	
217200	C3 CryoProbe 2.1 mm sharp for C3 CryoSystem	Diameter 2.1 mm, sharp tip Stem length 133 mm Handle length 140 mm Hose length 2 m	

Bibliography

- Hippocrates. Heracleitus on the Universe. Aphorisms, Vol4, Heinemann, London, 1931; 5:165, 7:201.
- Gruner OC. A Treatise On The Canon of Medicine of Avicenna. Luzac, London, 1930.
- Larrey D. Surgical Memoirs of the Campaigns of Russia, Germany and France. Carey & Lea, Philadelphia, 1832.
- Evans P. Cryoanalgesia. The application of low temperatures to nerves to produce anaesthesia or analgesia. Anaesthesia 1981; 36:1003-1013.
- Arnott, J. On Severe Cold or Congelation as a Remedy of Disease. London Medical Gazette, 1848.
- Arnott, J. On the treatment of cancer by the regulated application of an anesthetic temperature. J Churchill, London, 1851:32-54.
- Trendelenburg W. Über langdauernde Nervenanschaltung mit sicherer Regenerationsfähigkeit. Zeitschrift für die gesamte experimentelle Medizin 1917; 5:371-374.
- Amoils SP. The Joule Thomson Cryoprobe. Arch Ophthalmol 1967; 78:201-207.
- Lloyd JW, Barnard JDW, Glynn CJ. Cryoanalgesia. A new approach to pain relief. Lancet 1976; 2:932-934.
- Warfield CA. Cryoanalgesia: Freezing of peripheral nerves. Hospital Practice 1987:71-77.
- Barnard, D. The effects of extreme cold on sensory nerves. Ann R Coll Surg Engl 1980; 62:180-187.
- Evans PJ, Lloyd JW, Green CJ. Cryoanalgesia: the response to alteration of temperature and freeze cycle. Br J Anaesth 1981; 53:1121-1126.
- Evans PJ, Lloyd JW, Jack TM. Cryoanalgesia for intractable perineal pain. J R Soc Med 1981; 74:804-809.
- Myers RR, Powell HC, Heckman HM et al. Biophysical and pathological effects of cryogenic nerve lesion. Ann Neurol 1981; 10:478-485.
- Sunderland S. Nerves and Nerve Injuries. Edinburgh & London: Livingstone, 1968:180.
- Myers RR, Heckman HM, Powell HC. Axonal viability and the persistence of thermal hyperalgesia after partial freeze lesions of nerve. J Neurol Sci 1996; 139:28-38.

inomed Medizintechnik GmbH
Im Hausgruen 29
79312 Emmendingen (Germany)

Tel: +49 7641 9414-0
Fax : +49 7641 9414-94
info@inomed.com



→ [inomed.com](https://www.inomed.com)