

LG4



LG4 - a versatile radiofrequency generator for pain treatment



C3 - a dedicated cryoneurolysis system


Fire & Ice

inomed is the only company to offer both radiofrequency & cryoanalgesia systems



HEAT USED TO TREAT PAIN

The use of heat to treat pain has a long history that can be traced back to the early Egyptian civilisation over 3,000 years ago



Egyptians and other civilisations used heat to treat muscular spasms and pain

In 1866, Busch discovered that high temperatures cause tumours to shrink

In 1931, Kirschner described the use of radiofrequency for the treatment of trigeminal neuralgia

In 1960, Rees developed the technique of percutaneous rhizolysis

In 1980, Bogduk and Long developed percutaneous lumbar medial branch neurotomy

In 1998, Sluijter, Cosman, Rittman and van Kleef invented pulsed radiofrequency

2

In 1909, Nagelschmidt demonstrated that high-frequency energy causes tissue necrosis. He later introduced the first diathermy system

In 1950, Cosman and Arnoff developed the first commercial RF generator

In 1974, the first dorsal root radiofrequency was described

In 1975, Shealy described the first use of radiofrequency for lumbar facet pain

3


What is continuous radiofrequency (CRF) denervation and how does it work?

While heat has been used to treat pain for thousands of years, the use of radiofrequency to treat pain is a relatively new innovation. Radio waves are on the lower end of the electromagnetic spectrum and have the longest wavelengths and lowest energy of all types of radiation.

The term 'continuous radiofrequency' refers to the oscillation frequency of an alternating electric current; the frequency range used for radiofrequency denervation is between 250 and 500 kHz. The rapid oscillation of the radio waves causes ionic motion of the tissue directly adjacent to the active tip of the CRF electrode. The ionic motion leads to molecular friction and heating of the tissue within a predictable distance from the RF electrode. The heat causes coagulative necrosis of the nerve, which in turn results in a blockage of the nociceptive pathway.



In a monopolar CRF procedure, the patient is part of a closed loop circuit that includes the RF generator, electrode, and a large dispersing electrode.

Electrolytic effect caused by RF in human tissue





What is pulsed radiofrequency (PRF)?

Pulsed radiofrequency differs from continuous radiofrequency as it is a non-destructive treatment method where the temperature is kept below 42°, which is not hot enough to cause coagulative necrosis. Instead, the pulsed radiofrequency generates an electromagnetic field around the nerve that modulates the pain signal by altering the nerve's gene expression.




How can understanding the shape and size of the lesion improve the outcome?


The heat generated during a continuous radiofrequency procedure originates from the tissue rather than the electrode, which affects the distribution of heat within the tissue. The heat radiates laterally from the active tip of the electrode, but does not reach the area distal to the electrode tip.



This means that the electrode must be positioned parallel to the nerve, rather than perpendicular to it, in order to maximise the contact area with the nerve.



The heat lesion can be adjusted by altering certain parameters of the procedure.




In addition, a change in the size of the active tip also leads to an enlargement of the lesion.


	65 °C	70 °C	75 °C	80 °C
Diameter 22G (0.7 mm) Active tip 5 mm	1.4 mm	1.8 mm	2.1 mm	2.4 mm
Diameter 20G (0.9 mm) Active tip 5 mm	2.0 mm	2.3 mm	2.7 mm	3.1 mm

What type of injury is caused by the treatment and how does it differ from PRF and cryoanalgesia?

According to the Sunderland classification of nerve injury, a heat-induced lesion corresponds to a grade 3 injury.

A key difference between a lesion created by cryoanalgesia and one created by continuous radiofrequency is that the endoneurium is damaged during continuous radiofrequency, whereas it remains intact during cryoanalgesia. Damage to the endoneurium can lead to procedural neuritis and neuromas, which is not possible with lesions caused by cryoanalgesia. Pulsed radiofrequency is not included in the Sunderland classification, as the treatment process is based on neuromodulation rather than coagulative necrosis or Wallerian degeneration.





Indications of continuous radiofrequency (CRF) and pulsed radiofrequency (PRF)

CRF

- Intercostal
- Facet joint
- Genicular
- Gasserian ganglion
- Occipital neuralgia
- Chordotomy
- Sacroiliac joint

PRF

- Atypical facial pain
- Neuropathic pain in one or more divisions of the trigeminal nerve
- Postherpetic neuralgia
- Intractable cluster headache
- Dorsal root ganglion
- Hip joint



Although there are numerous overlaps in the indications, since pulsed radiofrequency does not use heat as part of the therapeutic process, it is mainly used in areas where motor and sensory functions need to be preserved. This also avoids some of the complications associated with thermal injuries, such as neuritis, weakness, or deafferentation pain.

What is the difference between radiofrequency and cryoanalgesia?

Continuous radiofrequency (CRF) achieves effective lesioning through the prolonged application of continuous electrical current, elevating tissue temperatures to 60–80 °C. This creates targeted neuroablative effects through coagulative modification of both cellular and acellular structures. Histologically, this results in precise changes within the axonal and connective tissue components of the nerve, including the endoneurium, perineurium, and epineurium. CRF is associated with reliable and sustained pain relief and is especially valued for its precision, which allows for focused treatment of pain pathways.

Pulsed radiofrequency (PRF) applies short bursts of radiofrequency energy that generate an electromagnetic field, influencing cellular activity in a way that positively modulates pain signal transmission. While the exact mechanisms are still under investigation, PRF is believed to affect gene expression, neuronal membrane function and cytokine profiles. PRF is non-destructive, preserving the structural integrity of the nerve and is particularly well-suited for situations requiring a gentler neuromodulatory approach. It offers flexibility and is well-tolerated, making it a favourable option in a variety of clinical scenarios.

Cryoanalgesia uses extreme cold to create a temporary and reversible block of pain signal transmission by forming an ice ball around the target nerve. This results in a controlled axonotmesis injury, producing a lesion that can be easily visualised via ultrasound. The technique creates a larger treatment area, allowing effective targeting of both pure sensory and mixed sensorimotor nerves without long-term motor impairment. Its versatility makes it highly valuable where precise and reversible nerve modulation is essential.

inomed Lesion Generator (LG) family

LG1

- 50 W generator
- Continuous radiofrequency: power- or temperature-controlled, or step mode
- Pulsed radiofrequency: temperature-, voltage-, or pulse width-controlled
- Sensory and motor stimulation
- 1 channel monopolar or 1 channel bipolar
- Continuous temperature monitoring
- Continuous impedance monitoring
- Contact quality monitoring
- Audio impedance
- Remote control

Art. No. 261000

LG2

- 50 W generator
- Continuous radiofrequency: power- or temperature-controlled, or step mode
- Pulsed radiofrequency: temperature-, voltage-, or pulse width-controlled
- Sensory and motor stimulation
- 2 channels monopolar, 2 channels monopolar as bipolar, 1 channel bipolar
- Continuous temperature monitoring
- Continuous impedance monitoring
- Audio impedance
- Contact quality monitoring
- Remote control

Art. No. 262000

LG2-sEEG

- 50 W generator (10 W for sEEG)
- 1 channel + sEEG
- Continuous temperature monitoring
- Continuous impedance monitoring
- Contact quality monitoring
- Audio impedance
- sEEG simulator test plug
- Remote control

Art. No. 263000

LG4

- 100 W generator
- Various monopolar and bipolar configuration options
- Sensory and motor stimulation
- Continuous radiofrequency: power- or temperature-controlled, or step mode
- Pulsed radiofrequency: temperature-, voltage-, or pulse width-controlled
- Continuous temperature monitoring
- Continuous impedance monitoring
- Contact quality monitoring
- Audio impedance
- Separate modes for different channels
- Remote control
- Foot switch control

Art. No. CR2020G



→ LG1



→ LG2



→ LG2-sEEG



→ LG4

Unique features of the LG4

- Large touchscreen display for intuitive control of the device
- Automatically controlled RF output with simultaneous or staggered starts
- Multiple lesion setups
- Neutral electrode contact quality monitoring to reduce the risk of skin burns
- Continuous impedance and temperature monitoring during all modes up to 4 channels
- Control via foot switch or remote control
- Patient database for treatment, detailed export of the protocol via USB stick
- Save up to 20 use scenarios per application















































Technology that works the way you do: pain accessories

Disposable cannulas / reusable (SuperLight) and disposable electrodes

Percutaneous RF cannulas

Straight und sharp cannulas

Art. No.	Description Type - diameter x active tip		
Length=50 mm			
240100	straight - 22G x 4 mm		
240109	straight - 22G x 7 mm		
Length=100 mm			
240101	straight - 22G x 2 mm		
240102	straight - 22G x 5 mm		
240114	straight - 22G x 7 mm		
240106	straight - 22G x 10 mm		
240111	straight - 20G x 2 mm		
240110	straight - 20G x 5 mm		
240112	straight - 20G x 10 mm		
240120	straight - 17G x 2 mm		
240121	straight - 17G x 5 mm		
240122	straight - 17G x 7 mm		
240123	straight - 17G x 10 mm		
Length=115 mm			
240140	disposable cannula - 18G x 2 mm		
240141	disposable cannula - 18G x 5 mm		
240142	disposable cannula - 18G x 7 mm		
240143	disposable cannula - 18G x 10 mm		
Length=150 mm			
240103	straight - 20G x 2 mm		
240104	straight - 20G x 5 mm		
240105	straight - 20G x 7 mm		
240108	straight - 20G x 10 mm		
240107	straight - 20G x 15 mm		

Fitting electrodes for cannulas

RF pain electrodes

Art. No.	Description
Length=50 mm	
260007	disposable
260006	reusable (SuperLight)
Length=100 mm	
260012	disposable
260011	reusable (SuperLight)
Length=115 mm	
260126	disposable
260121	reusable (SuperLight)
Length=150 mm	
260017	disposable
260016	reusable (SuperLight)



Bibliography

1. Wray J, Dixon B, Przkora R. Radiofrequency ablation. In: StatPearls. Treasure Island (FL): National Library of Medicine; 2023 Jun 12.
2. Patel P, Danish S. Radiofrequency ablation techniques. ScienceDirect.
3. Radiofrequency ablation: mechanism of action and devices. J Vasc Interv Radiol. 2010 Aug;21(8):179-189.
4. Mechanism of action of radiofrequency ablation. Springer Nature.
5. Essentials of radiofrequency ablation of the spine and joints. Springer Nature.
6. Kirschner M. Zur Kryochirurgie. Arch Klin Chir. 1931;167:761-768.
7. A history of the development of radiofrequency neurotomy. J Pain Res.
8. Choi EJ, Choi YM, Jang EJ, Kim JY, Kim TK, Kim KH. Neural ablation and regeneration in pain practice. Korean J Pain. 2016 Jan;29(1):3-11.
9. McCormick Z, Hurley R. The evolution of radiofrequency denervation for pain indications. Pain Med. 2021.
10. Cosman ER Sr, Cosman ER Jr. Electric and thermal field effects in tissue around radiofrequency electrodes. Pain Med. 2005.
11. Erdine S, Billir A, Cosman ER Sr, Cosman ER Jr. Ultrastructural changes in axons following exposure to pulsed radiofrequency fields. Pain Pract. 2009;9(6):407-417.
12. Cosman ER Jr, Dolensky JR, Hoffman RA. Factors that affect radiofrequency heat lesion size. Pain Med. 2014.
13. Sunderland S. Nerves and nerve injuries. Edinburgh: Livingstone; 1968.

inomed Medizintechnik GmbH
Im Hausgruen 29
79312 Emmendingen
(Germany)

Phone: +49 7641 9414-0
Fax: +49 7641 9414-94
info@inomed.com



→ [inomed.com](https://www.inomed.com)

D030307-EN / Content subject to change © copyright by inomed / 05.2026 / REV 1.0