

star

RF Electrode_Fixed

inomed 


ADVANCED TECHNOLOGY IN RFA



- Internally cooled electrode
- Wide selection of needle gauges, lengths and active tips
- Lightweight handle for easy manoeuvring
- The first thyroid-dedicated electrode on the market

 **STAR**med

CE
1639



Specification of star RF Electrode - Fixed

Gauge	Length	Active Tip
15 G	150 mm	15 mm, 20 mm, 25 mm, 30 mm, 40 mm
	200 mm	20 mm, 25 mm, 30 mm, 40 mm
	250 mm, 350 mm	20 mm, 30 mm, 40 mm
17 G	150 mm	7 mm, 10 mm, 15 mm, 20 mm, 25 mm, 30 mm, 40 mm
	200 mm	10 mm, 15 mm, 20 mm, 25 mm, 30 mm, 40 mm
	250 mm, 350 mm	10 mm, 20 mm, 30 mm, 40 mm
18 G	70 mm	4 mm, 5 mm, 7 mm, 10 mm, 15 mm
	100 mm	5 mm, 7 mm, 10 mm, 15 mm, 20 mm
	150 mm	5 mm, 7 mm, 10 mm, 15 mm, 20 mm, 30 mm
19 G	70 mm, 100 mm	4 mm, 5 mm, 7 mm, 10 mm, 15 mm, 20 mm

Manufacturer



STARmed Co., Ltd.
#B-14th Floor, Daebang Triplaon
158 Haneulmaeul-ro, Ilsandong-gu
Goyang-si Gyeonggi-do
10355 Korea

Distributor



inomed Medizintechnik GmbH
Im Hausgruen 29
79312 Emmendingen (GERMANY)
Tel. +49 7641 9414-0 · Fax +49 7641 9414-94
info@inomed.com · www.inomed.com