inomed

C2 Xplore





One step ahead –

digital transformation in IONM

Simplification

by digitalization and automation



C2Xplore



The ultimate in user-friendliness

Automated and guided applications
with an intuitive touch function:
gesture-controlled comment function
and step-by-step guidance for clarity

Maximum user-friendliness
with intuitive and clear display

Digital health record:
simple and secure documentation
during surgery supported by
simplified data export
and an integrated database

Flexibility for cost-effectiveness

1 device › 2 application fields



Thyroid surgery

To locate and monitor the recurrent laryngeal nerve and vagus nerve



ENT surgery

To locate and monitor the facial nerve

C2Xplore

Cost efficiency by versatility


Software flexibility
with up to 2 different applications in 1 device

Ease-of-use
a range of applications – yet just one C2 Xplore

New hardware
with a sustainable concept, proven and tested for more than 10 years



Native OR integration



* The installation is project-based and needs to be evaluated within your hospital IT department.

C2Xplore

Maximum interconnection

Facilitated data export

Digital integration in the OR
Data export by Wi-Fi to the clinic network by HL7*
PDF export or PostScript printer

Presence in the OR
Connection to various screens or wireless with a tablet

inoView V2 connection to any tablet
(from Android 8.0 and iOS 12.1)
to suit all personal preferences

Control comments and settings
inside the sterile area with inoView V2

2D LED scanner
for super-fast data acquisition of patient data

Reports printed immediately
Direct printer or network printer

Safety not only for the nerves
No data loss thanks to the large internal storage



* The installation is project-based and needs to be evaluated within your hospital IT department.



Connectivity – Inside and outside the OR

Amplify your research and networking

Join the international network of inomed users to share work and cooperation

Evaluate surgery with the C2 Xplore's dedicated analysis function



C2 Xplore



Support at your fingertips

More confidence through immediate support

Remote support for instant help and diagnosis in the OR for more than 10 years

International IONM exchange



inomed

Broaden your knowledge and connect with international experts



Online webinars and training with other users to broaden IONM knowledge, and OR visits for demonstration of cutting-edge technology. Broad mix of offline and online education for maximum efficacy

C2Xplore


inomed

>>> We manufacture products which are used in more than 5,000 hospitals in 100 countries to perform around 500,000 operations per year.



A single provider for all IOM needs

We offer complete intraoperative neuromonitoring solutions. Our extensive range of accessories combined with direct, personal contact and fast, on-site service anywhere in the world guarantee that you receive support you can rely on. Our courses and workshops ensure a consistently high standard of service and quality.



Continuous innovation and new fields of use

As a technology leader, we work with renowned experts worldwide on an extensive program of research and development to continuously create new fields of use and product innovations.

Developed and manufactured to German quality standards

With responsibility, attention to detail and a pioneering spirit, since 1991 we manufacture products for the benefit of patients and doctors.

D030116-US / EN / Content subject to change © copyright by inomed / 01.2023 / REV1.0

inomed

inomed Inc.
1 Westbrook Corporate Center
Suite 830
Westchester, IL 60154
USA

Tel.: +1 708 607 3250
Fax: +1 708 607 3259
info.us@inomed.com
www.us.inomed.com